

# NOBILTRAY

## — Description

NOBILTRAY is a preformed light curing tray indicated for making individual tray and can also work as temporary base plate in dental labs. As light curing material, it can be fast cured in 3 minutes in a unit. The product is characterized by easy operation, high precision, low deformation rate, high plasticity and pleasant mint flavour.



## — Application

- Making individual tray



| Types / Features | Light Curing Tray  | Self-cured resin Tray   | Metal Tray  |
|------------------|--|---|---|
| Operation        | preformed → ready for use and time-saving                                    | need mixing and kneading first → complicated and time-consuming                         | -   |
| Plasticity       | soft but with certain strength → easy to trim and shape when opening windows | too soft to work with in dough stage → difficult to trim and shape when opening windows | hard and rigid metal material → difficult to open window                      |
| Accuracy         | high hardness after curing → low deformation rate ensuring accuracy          | existence of contraction → unstable and inaccurate                                      | specific and uniform size → fail to match different oral forms and inaccurate |
| Smell            | mint flavour → pleasant smell  | residual stimulative monomer → pungent smell  | -   |

- Making temporary base plate

| Types / Features | Light Curing Base Plate   | Wax base Plate                                      |
|------------------|---|---|
| Operation        | preformed → ready for use and time-saving                           | need heating first → time-consuming                 |
| Accuracy         | high hardness after curing → low deformation rate ensuring accuracy | vulnerable to temperature → unstable and inaccurate |



## — Main Advantages

- Preformed → easy handling
- 3-minute fast curing → quick curing
- High hardness and fit different oral forms → accurate
- Soft → high plasticity
- Mint flavour → pleasant smell
- Various choices → pink or blue, thicker or thinner, all available

## — Specification

Standard package



50pcs / box

Sample package



2pcs / bag



Follow @VinciSmile

VINCISMILE GROUP LLC

Add: 20524 Carrey Rd. Walnut, CA 91789, USA

Email: marketing@vincismile.com

Web: www.vincismile.com